CASTCO HARD**Co-Cr Dental Ceramic Alloy Instruction --- Nikelsiz Kobalt Krom Seramik Metali Kullanma Talimatı**

TL.09.33 / Y.T.:06.07.2011 / R.T.:06.05.2014 / Rev.No:02

TS EN ISO 22674 Standartlarına uygundur

Chrome Cobalt Ceramic Alloy – Nickel Free, Beryllium Free**HARD**

T e c h n i c a l D a t a	
Melting Range / Ergime Derecesi	1350-1400 °C
DuctileYield (A5) / Yield Modülü (A5)	% 6
Modülüs to Elasticity / Esneklik Ölçüsü	210,000 (MPa)
Elongation / Elangasyon	390 (MPa)
Density / Yoğunluk	8.5 (g/cm3)
Vickers Hardness / Vickers Sertliği	360
Coefficient of the liner Expansion / Liner Genleşme Kat Sayısı	14-14.2 x 10-6 K-1
Tensile strength [RM]/ Çekme (Gerilme)	720
C o m p o s i t i o n	
Cobalt /Kobalt	62.5 %
Chrome / Krom	24 %
Niobium /Niobyum	1.5 %
Wolfram / Wolfram	8.5 %
Vanadium / Monibden /Other	V 1.31 % / Mo 1.53 % / Si 0.64% Fe 0.82 % C 0.019 % Mn 0.031 %

Instruction:

Modelling: Minimum preparation for metal and Porcelain Single crowns : 1.3 – 1.5 mm Bridgework : 2.5 – 2.0 mm

Before the wax-up the preparations should be coated with die varnish. The size of the metal crown should be reduced to correspond to that of the final prosthetic crown, missing portions of the teeth. A cross-sectional wall thickness of 0.4 mm at the accusal surface will ensure a successful casting. Uniform porcelain thickness ensures stress-free bonding. In order to avoid chipping of the porcelain, the lingual metal margin should be placed below the incisal region. Contact surfaces provided for soldering and welding. Rounding of the metal-porcelain interface ensures cleaner margins and an esthetic colour. The outstanding properties permit very thin connections to be used between the crowns and pontics. During later firing process the individual teeth can be seated clear of the adjacent teeth. This avoids a block appearance. Sprue System: Minimum metal thickness (after grinding) : for ceramic veneering 0.3 mm for acrylic veneering with retention pearls 0.3 mm. Shape groove in cervical and palatal area. Avoid sharp edges. In the case of work on long bridges, divide plastic hollow stick between the front teeth and in the area around the canine teeth. Investing / Preheating : Use phosphate-bonded crown and bridges investment materials. (DENTindex Investment) Follow processing instructions. Preheating temperature 900-1000 °C Melting / Casting: Do not overheat alloy. Use only clean crucibles, one crucible per alloy. It is to enable an exact identification of each charge cast new metal only. In case of re-casting: only re-cast identical alloys. Blast old material. Add at least 50 % of new material. Use only ceramic crucibles. Moment for Casting: Vacuum pressure casting with induction heating and centrifugal casting with induction heating. When the last solid component has submerged completely in the melt, continue to heat for 0 to 12 seconds depending on the induction capacity of the casting unit, then trigger casting. Follow operating instructions for induction furnaces. Flame centrifugal casting: Cast when the last solid component has submerged completely in the melt and the melt clearly moves due to the flame pressure. After deflasking: Blast residual investment material. Use fine carbide, ceramically bonded stones or sintered diamond milling tools for finishing. Final Work : Blast visible metal surfaces, then blast-polish external surfaces. Rubber polish with DENTindex rubber polishers. Clean thoroughly steam clean or boil in aqua dest. Secondary Effects: Such as allergies to contents of the alloy or electrochemically based reactions may very rarely occur. Reactions: In case of known incompatibilities and allergies to contents of the alloy. Warranty : Whether given verbally, in writing or by practical instructions, our recommendations for use are based upon our own experience and trials and can only be considered as standard values. Our products are subject to a constant further development. Therefore alterations in construction and composition are reserved.

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Marking: According to your request**Packing : Leather Pocket, Plastic Pocket, Wood Pocket****Diameter: 8 mm x 17 mm****Code: 11193-C**

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UĞURLU DIŞÇILIK GIDA ENDÜSTRİ & İNŞAAT ÜRÜNLERİ SAN.TİC.LTD.ŞTİ